

Los Angeles County Public Health Laboratory 12750 Erickson Avenue Downey, CA 90242 Phone (562) 658-1300 Fax (562) 401-5999

**Arbovirus Real-Time PCR (8/5/16)** 

|                                   | Albovilus Real-Tille FCR (8/3/10)   |
|-----------------------------------|---|
| Other Name(s)                     | CDC Trioplex RT-PCR (Dengue, Chikungunya, Zika virus)   |
| Sunquest Test Code                | TRPCR   |
| Pre-Approval Required             | No  |
| Supplemental Information          | Los Angeles County PHL Test Request Form  |
| and Required Form(s)              | http://www.publichealth.lacounty.gov/lab/labforms.htm   |
| Required Specimen                 | Serum collected within 7 days of symptom onset  |
| Type(s)                           | Urine (random) collected within 14 days of symptom onset  |
|                                   | CSF, amniotic fluid, tissue, cord blood by consultation only; CSF, urine and amniotic fluid each must be tested collected alongside a patient-matched serum specimen.   |
| Minimum Volume                    | 707   |
| Required                          | PCR can detect Dengue, Chikungunya, Zika virus  |
|                                   | Serum 5-7 mL<br>CSF 2mL   |
|                                   | CSF ZIIIL   |
|                                   | PCR can detect Zika only  |
|                                   | Urine 10-20 mL  |
|                                   | Amniotic fluid 5-10 mL  |
|                                   | Cord blood 5-7 mL   |
|                                   | Tissue 1cm <sup>3</sup>   |
| Storage/Transport                 | Store specimens at 4-8°C and ship to laboratory immediately on cold pack.   |
| Conditions                        | Ship formalin fixed samples at room temperature.  |
| Transport Medium                  | Transport medium not applicable.  |
| Specimen Labeling                 | Test subject to CLIA regulations and requires two unique patient identifiers on the primary   |
|                                   | specimen container and the test requisition including patient ID, date of collection, submitter information, and specimen ID number. The identifiers must be clearly labeled on specimen and must match to information on the requisition form. |
| Shipping Instructions and         | Contact public health laboratory for assistance with courier at 562-658-1460 during   |
| Specimen Handling<br>Requirements | business hours. After hours, contact the County Operator at 213-974-1234 and ask for the Public Health Laboratories Director.   |
| Test Methodology                  | Real-time PCR   |
| Turnaround Time                   | 7 business days for urine, csf, serum/cord blood, and amniotic fluid Tissue specimen types are sent out testing to CDC  |
| Interferences & Limitations       | Negative PCR results do not necessarily rule out Zika virus. Additional specimen may be required for IgM serology.  |
| Additional Information            | Emergency Use Authorization Test  |
| Reference Range                   | Not detected  |
| CPT Code(s)                       | 87801   |