

FENTANYL TEST STRIPS CAN SAVE YOUR LIFE!



What are fentanyl test strips?

- Test strips can tell you if drugs contain fentanyl.
- They do not tell you how much fentanyl there is, or how strong it is. But they can help prevent overdose if used correctly.

How do they work?

- Fentanyl test strips are dipped into drug residue dissolved in water. Within minutes, a person can know whether the drug contains fentanyl.

Where can I get test strips?

- Visit your local health center to get **FREE** test strips.
- You can also visit ph.lacounty.gov/clinics for more information.

HOW TO USE FENTANYL TEST STRIPS

Step 1: Add water and mix


- For Meth, MDMA and ecstasy, use 1 teaspoon of water for every 10mg of crystal or powder you are testing.
- For all other drugs, use a half teaspoon of water.

Step 2: Use the test strip

- Place the test strip with the wavy side down in the water. Let the strip absorb the water for 15 seconds.
- Take the strip out of the water and place it on a flat surface for two minutes.

Step 3: See results

- If fentanyl is found in the drug, one line will be present.
- If fentanyl is not found, two lines will be present.



Learn More:



SCAN ME