



COUNTY OF LOS ANGELES  
DEPARTMENT OF PUBLIC HEALTH

PUBLIC HEALTH LABORATORY  
12750 ERICKSON AVENUE  
DOWNEY, CA 90242

PHONE (562) 658-1300  
FAX (562) 401-5999

PLACE BARCODE LABEL HERE

CT/GC NAAT TEST REQUEST FORM

FOR RYAN WHITE HIV CLINICS FUNDED BY  
THE DIVISION OF HIV AND STD PROGRAMS

PATIENT NAME (LAST, FIRST)				DATE COLLECTED		TIME COLLECTED	
PATIENT ADDRESS				SUBMITTER / CLINIC			
PATIENT / PARTICIPANT MRN#		M <input type="checkbox"/>	F <input type="checkbox"/>	OTHER <input type="checkbox"/>	DATE OF BIRTH		PROGRAM ACCOUNT <b>TMP</b>
REQUESTING PHYSICIAN		PHYSICIAN #		SPECIMEN SOURCE		DATE / TIME RECEIVED (FOR LAB USE ONLY)	

THERAPEUTIC MONITORING PROGRAM ELIGIBILITY CERTIFICATION

THIS CLIENT HAS AN ANNUAL ADJUSTED GROSS INCOME OF LESS THAN \$50,000 AND DOES NOT HAVE MEDI-CAL, MEDICARE, OR OTHER THIRD PARTY INSURANCE. CITY/COUNTY INDIGENT PROGRAMS DO NOT FALL UNDER THE CATEGORY OF THIRD PARTY INSURANCE.

YES INITIALS: \_\_\_\_\_ CHECK COVERAGE TYPE:  RYAN WHITE  HEALTHY WAY L.A. (UNTIL 12/31/13)

TEST	SPECIMEN TYPE / REQUIREMENTS	SPECIMEN STORAGE / TRANSPORT
<input type="checkbox"/> C. TRACHOMATIS/N.GONORRHOEAE (CT/GC) NAAT	<input type="checkbox"/> URINE- COLLECT A MINIMUM OF 2ML URINE WITH APTIMA URINE SPECIMEN COLLECTION KIT (YELLOW) <input type="checkbox"/> RECTAL- COLLECT WITH APTIMA VAGINAL SWAB SPECIMEN COLLECTION KIT (ORANGE) <input type="checkbox"/> VAGINAL- COLLECT WITH APTIMA VAGINAL SWAB SPECIMEN COLLECTION KIT (ORANGE)	STORE COLLECTION KITS AT 15-30°C STORE AND TRANSPORT SPECIMENS AT 2-30°C SEND SPECIMENS TO LABORATORY AS SOON AS POSSIBLE FOR TESTING REJECTION CRITERIA INCLUDES SPECIMEN WITH INCORRECT/INADEQUATE PATIENT IDENTIFICATION; UNLABELED SPECIMEN; EXPIRED COLLECTION KIT; INAPPROPRIATE TRANSPORT CONDITIONS; LEAKING SPECIMEN; URINE SPECIMEN RECEIVED >30 DAYS AFTER COLLECTION; SWAB SPECIMEN RECEIVED >60 DAYS AFTER COLLECTION; WHITE CLEANING SWAB IN SPECIMEN COLLECTION TUBE; SWAB SPECIMEN SUBMITTED WITHOUT SWAB; MULTIPLE SWABS IN SPECIMEN COLLECTION TUBE; INCORRECT SPECIMEN VOLUME; INAPPROPRIATE COLLECTION KIT FOR SPECIMEN SOURCE
<input type="checkbox"/> N.GONORRHOEAE (GC ONLY) NAAT	<input type="checkbox"/> THROAT- COLLECT WITH APTIMA VAGINAL SWAB SPECIMEN COLLECTION KIT (ORANGE)	